
Specification

Article	Specification			
Strength	1~6 Bar, 0.1 Bar step			
Frequency	1~21Hz			
Treatment Mode	M1~M20			
Treatment Shot	100~5000, continuous, 100 step			
Channel	1 or 2Ch(option)			
Applicator	6, 10, 15, 20, 35mm(option)			
Interface	8" LCD, Touch Screen, Jog Shuttle			
Dimension	467 x 501 x 993mm			





Factory

Seoul Branch 4F,224,Baekjegobun-ro,Songpa-gu,Korea

Tel:+82-1588-7395 Fax:+82-2-418-0986


#301~303,187,Techno 2ro,Yuseong-gu,Daejeon-si,Korea

Tel:+82-42-934-5560~1 Fax:+82-42-934-5562

s a l u s - R/S///W

The World's Best Digital Control Method
Radial Shock Wave Therapy

Salus-RSWT



1. High Power(1~6 Bar), Differentiated Frequency (1~21Hz) RSWT

6mm, 10mm, 15mm, 20mm, 35mm

6Bar 21Hz, 2Channel Radial Shock Wave

Salus-RSWT provides various applicators.

- 2. 8" LCD User Interface (Simple operation, Various Treatment Mode)
- 3. Fine-Tune Output(0.1 Bar)
- 4. Low Noise Pneumatic Compression Device
- 5. Low Maintenance Cost
- 6. Luxurious Design
- 7. High Durable Solenoid Valve
- 8. 20 Memory Operation Mode (Name Tag Function)

Various Applicators

6mm	10mm	15mm	20mm	35mm(option)



" Authentic Radial Shock Wave

Radial Shock Wave Feature

